



Product Main Number 1000021307

ABENA®

Shaped pad, ABENA San Premium, 5, Purple colour code

- ✓ High performance core with Top-Dry Technology
- ✓ 360° leakage protection preventing leakage all around
- ✓ Wetness indicator with graduation scale
- ✓ Support healthy skin balance
- ✓ Cost-effective solution
- ✓ 100 % breathability for healthy skin and optimum comfort
- ✓ Anatomically shaped for better fit
- ✓ Odour Control



Product Description

ABENA San 5-12 is our premium range of anatomically shaped pads for moderate to heavy incontinence. The high performance core with Top-Dry Technology supports optimal skin care by keeping the exposed area comfortable and dry even after multiple voids. The fully breathable pads are fitted with standing leakage barriers combined with front and rear built-in cross barriers offering 360° leakage protection. When combined, these functionalities help protect the skin and maintain a natural skin balance. The Odour Control system minimises unpleasant smells. The wetness indicator with graduation scale guides the user and carer to the product with the right absorption level and product changes. The dermatologically tested products have been given the Nordic Swan Ecolabel. ABENA San in sizes 5-12 should always be used with special fixation pants, like our fixation pants, to ensure a tight fit and reduce the risk of leakage. ABENA's two-piece systems, consisting of fixation pants and an incontinence pad, are a cost-effective solution.



Specifications

Base name	Shaped pad
Brand	ABENA
Sub-Brand	San Premium
Size	5
Colour	White
Features	Purple colour code
Material	ECF fluff, PE, breathable textile backsheet (BTBS), elastane, nonwoven, PP, resin, SAP
Length/depth	54 cm
Width	28 cm
Tolerance	+/- 10%
Absorption Rothwell	1200 ml
Gender	Unisex
Certificates	CE. Nordic Swan Ecolabel. MD. UKCA.
Nordic Swan Ecolabel License-No.	5023 0001
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745
Safety Instructions And Warnings	Keep out of reach of children.
Shelf Life	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations.

Packaging Disposal Instructions

Can be recycled or incinerated.

Instructions for use/application

Bladder protection pads for moderate to heavy incontinence. Easy product size recognition by colour codes, names and numbering systems. ABENA San products are suitable for all needs and body shapes as they are worn between the legs and not around the stomach and hip. ABENA San 5-12 should always be used with special fixation undergarments, such as our fixation pants. The fixation pants ensure a tight fit and minimise the risk of leakage. Pull the fixation pants half-way up the legs. Fold the pad lengthways to create a bowl shape. Pass the pad between the legs from front to back. Position the upper rear end of the pad at the tailbone. Fan out the back of the pad horizontally to create a channel shape between the legs. Use the wetness indicator to check if the product is applied accurately in the centre. Fan out the pad horizontally at the front. Pull up the fixation pants into position.

Packaging data

Unit	Contains	Length	Width	Height	EAN
car	4 pck	40.1 cm	40 cm	37.7 cm	5713571000571
pck	36 pcs	36 cm	16 cm	27 cm	5713571000564



The Swan label, also known as the Nordic Swan Ecolabel, is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material production, use, documentation and safety.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products and medical equipment.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



Leakage barriers in the leg cut and around the absorbent core help prevent leakages and provide an optimal fit.



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



The unique structure of our high performance core is designed for fast dispersion and encapsulation of urine, which results in a dry surface even after multiple voids.



Odor control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal Odor arises from the product.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



Wetness indicator is a liquid-sensitive strip that changes colour or fades out indicating when it is time to change the product.



A combination of 100% breathability, Top-Dry and leakage protection provides a dry surface and healthy skin even after multiple voids maintaining the skin's natural pH balance.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 12 drops.



The product does not contain latex.



Unique Top-Dry Technology that provides a rapid absorption and ensures a dry surface even after multiple voids.



The bag with this symbol is made with at least 60% PCR-plastic materials (post-consumer recycled), previously used by consumers. This allows us to reduce our plastic waste and to help take care of the environment.

