



Product Main Number 1000021314

ABENA®

shaped pad, ABENA San Premium Special, F4, Light green colour code

- ✓ For faecal incontinence
- ✓ Effective core technology with Top-Dry technology
- ✓ Unique leakage barriers with built-in pockets
- ✓ Supports healthy skin
- ✓ Odour Control minimises unpleasant smells
- ✓ 360-degree leakage barriers that prevent leakage all around
- ✓ Cost-effective solution
- ✓ 100% breathability for healthy skin and optimal comfort



Product Description

ABENA San Special is our anatomically shaped pad designed specifically for faecal incontinence. The pad can be used for both faecal and urinary incontinence and has a textile-like light-green backing foil for optimal comfort and discretion. The product has extra high leakage barriers that ensure that faeces stay in the nappy. The effective core with Top-Dry technology provides fast absorption and a dry surface, even after repeated urination, ensuring optimal skin balance and increasing comfort. The fully breathable pads are equipped with standing leakage barriers combined with built-in cross barriers at the front and back, providing 360-degree leakage protection. The combination of these functionalities helps to protect the skin and maintain a natural skin balance. The odour system minimises the risk of unpleasant smells. The wetness indicator helps the user and caregiver to find the product with the right absorption level and the right time to change the product. The dermatologically tested products carry the Nordic Swan Ecolabel. ABENA San Special should always be used with special fixation underwear, such as our ABENA Fix products. The fixation pants ensure a snug fit and reduce the risk of leaks.



Specifications

Base name	Shaped pad
Brand	ABENA
Sub-Brand	San Premium Special
Size	F4
Colour	White
Features	Light green colour code
Material	ECF fluff, PE, breathable textile backsheet (BTBS), elastane, nonwoven, PP, resin, SAP
Length/depth	73 cm
Width	37 cm
Tolerance	+/- 10%
Absorption Rothwell	2000 ml
Gender	Unisex
Certificates	CE. MD. Nordic Swan Ecolabel. UKCA.
Nordic Swan Ecolabel License-No.	5023 0001
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745
Safety Instructions And Warnings	Keep out of reach of children.
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations.

Packaging Disposal Instructions

Can be recycled or incinerated.

Instructions for use/application

Anatomically shaped diaper designed for combined faecal and urinary incontinence. Easy product recognition using coloured backing foil and printed name. ABENA San Special should always be used with special fixation underwear, such as our ABENA Fix products. The fixation pants ensure that the pad is kept close to the body and reduces the risk of leakage. Pull the fixation pants up to the middle of your legs. Fold the pad lengthwise to create a bowl shape with the back facing outwards. Pass the pad between your legs from front to back. Place the top edge of the pad at the tailbone. Spread the back of the diaper horizontally to create a cup shape between the legs. Use the wet indicator to check if the product is positioned exactly in the centre. Fold the pad out horizontally in front. Pull the fixation pants up into place. See our fitting instructions for detailed information.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cil	120 pcs	60 cm	40 cm	45.2 cm	5713571000717
pck	30 pcs	34.5 cm	20 cm	29 cm	5713571000700
pcs	1 pcs				



The Swan label, also known as the Nordic Swan Ecolabel, is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material production, use, documentation and safety.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products and medical equipment.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



Leakage barriers in the leg cut and around the absorbent core help prevent leakages and provide an optimal fit.



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



The unique structure of our high performance core is designed for fast dispersion and encapsulation of urine, which results in a dry surface even after multiple voids.



Odor control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal Odor arises from the product.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



Wetness indicator is a liquid-sensitive strip that changes colour or fades out indicating when it is time to change the product.



A combination of 100% breathability, Top-Dry and leakage protection provides a dry surface and healthy skin even after multiple voids maintaining the skin's natural pH balance.



The bag with this symbol is made with at least 60% PCR-plastic materials (post-consumer recycled), previously used by consumers. This allows us to reduce our plastic waste and to help take care of the environment.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 12 drops.



The product does not contain latex.



Unique Top-Dry Technology that provides a rapid absorption and ensures a dry surface even after multiple voids.

