



Product Main Number 1999915572



Protective sheet, ABENA Protective Sheet, 2 layers, 170x80cm

- ✓ High absorbency
- ✓ Peace of mind and protection around the clock
- ✓ Soft and comfortable material



Product Description

ABENA Protective Sheet is used primarily for bed and couch protection to ensure user peace of mind in any situation. The single use sheet has an absorbent surface and a liquid safe back sheet, which, combined with the sealed sides, allow the sheeting to provide optimum leakage protection. In addition to helping provide peace of mind when sitting or lying down, the use of protective sheets also helps reduce the risk of contamination. At the same time, the sheet ensures increased hygiene. ABENA Protective Sheets are available in different sizes.



Specifications

| | |
|----------------------------------|---|
| Base name | Protective sheet |
| Brand | ABENA |
| Sub-Brand | Protective Sheet |
| Secondary color | Light blue |
| Single or multiple use | Single use |
| Material | PE, paper |
| Length/depth | 170 cm |
| Width | 80 cm |
| Tolerance | +/- 3 cm |
| Layer quantity | 2 layers |
| Absorption Rothwell | 250 ml |
| Certificates | CE. MD. UKCA. |
| CE Class (Medical Devices) | Class I |
| Product or test standards | ISO 11948-1 |
| Directives, regulations and acts | MDR (EU) 2017/745 |
| Shelf Life | 5 years |
| Storage Instructions | Store dry, room temperature and no direct sunlight. |
| Product Disposal Instructions | Can be disposed of with normal household waste sorted according to local regulations. If contamination is suspected, the product must be discarded as clinical waste. |
| Packaging Disposal Instructions | Can be recycled or incinerated. |

Instructions for use/application

Ensure that the protective sheet is placed on the bed or chair with the back sheet down. Check that the sheet is fully spread out and smoothed to remove wrinkles. The sheet can be disposed of with normal household waste.

Packaging data

| Unit | Contains | Length | Width | Height | EAN |
|------|----------|---------|---------|--------|---------------|
| car | 4 pck | 45.9 cm | 23.6 cm | 20 cm | 5713571007204 |
| pck | 25 pcs | 47.2 cm | 26.4 cm | 7 cm | 5713571007198 |



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products and medical equipment.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



The product is only for single-use.



Barriers lengthwise of the product can help reduce incontinence related leakages onto beds and other furniture or surfaces.

