



Product Main Number 4103

ABENA ABRI-FIX Leaf


### Fixation pants, ABENA Abri-Fix Leaf, L, white, with legs, brown colour code

- ✓ **Comfort:** breathable soft polyamide material with legs
- ✓ **No seams** preventing pressure marks and skin irritation
- ✓ **Leakage security** provided by elastic knitting construction and fixation zones
- ✓ **Underwear-like look** for discretion
- ✓ **Easy product identification** with colour coding, washing instruction and material composition
- ✓ **Cost-in-use:** high cost effectiveness through a longer usage
- ✓ **Broad range** (XS-XXL)



### Product Description

Abri-Fix Leaf is an elastic fixation pant with legs made of breathable, soft polyamide material featuring stabilising zones. The pants look like normal underwear. The elastic material and knitted construction holds both small and large pads tight to the body preventing leaks. The seamless knitting technology – an innovative concept in producing clothing almost free of seams - prevents pressure marks and skin irritation. The products can be washed 100 times, therefore have a high cost effectiveness through a longer usage. Abri-Fix Leaf has coloured text representing size, washing instructions and material composition sewn in the waistband. 95% polyamide, 5% elastane



### Specifications

Base name	Fixation pants
Brand	ABENA
Sub-Brand	Abri-Fix Leaf
Size	L
Colour	White
Properties	With legs
Features	Brown colour code
Ingredients / Composition	95% polyamide, 5% elastane
Length/depth	330 mm
Width	110 mm
Waist circumference	85-115 cm
Thigh circumference	55-75 cm
Hip size	85-135 cm
Weight, net	25 g
Gender	Unisex
Certificates	CE. OEKO-TEX. MD.
CE Class (Medical Devices)	Class I
Directives, regulations and acts	MDR (EU) 2017/745
Shelf Life	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Washing Instructions	Wash at 60°C. Tumble dry on gentle program. Washable app. 100 times. The pants should be washed in laundry bag, which helps to extend shelf life. Use of bleaching agents may compromise product shelf life. Shrinkage approx. 7%.

<b>Product Disposal Instructions</b>	Can be disposed of with normal household waste sorted according to local regulations.
<b>Packaging Disposal Instructions</b>	Can be recycled or incinerated.

### Instructions for use/application

Abri-Fix Leaf should be worn as part of the 2 piece absorbent incontinence solution, for fixation of an incontinence pad. Fixation pants are wider at the crotch area and more elastic than normal underwear. This way they are making everyday-life easier by providing high comfort and security.

### Packaging data

Unit	Contains	Length	Width	Height	EAN
cil	30 pcs	34 cm	17 cm	23.8 cm	5703538124150
pck	3 pcs	32 cm	20 cm	3 cm	5703538124136
pcs	1 pcs				



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products and medical equipment.



OEKO-TEX, also known as "Confidence in Textiles", is a health and environmental label used throughout the world. The certificate ensures that textiles do not contain substances that are harmful to health or the environment.



The product does not contain latex.



The product can be washed at 60 degrees or below.



Do not bleach the product.



The product can be dried in the tumble dryer at low temperatures.



Do not iron the product.



Do not dry clean the product.



The product can be washed up to 100 times.

