



Product Main Number 491910



## Wash Glove, Molton, ABENA, 23x16cm, white, ultrasonic sealed, single use

- ✓ The soft material cares for and protects delicate and sensitive skin
- ✓ Extra high absorbency and good strength when wet
- ✓ Practical and hygienic solution
- ✓ Suitable for full-body washing
- ✓ Ideal for washing without water

### Product Description

ABENA's Wash Glove for personal hygiene is made of strong material and designed for single-use. The glove's embossed surface makes it easy to remove impurities and bacteria from the skin. The glove is intended for a full cleansing, and its soft and comfortable material makes it particularly suitable for washing exposed areas and delicate skin. The glove is strong and absorbent, even when wet. The practical opening in the wash glove makes it easy to put on. We also recommend using an ABENA nitrile glove under the wash glove for extra protection against bacteria and viruses. Our wash glove is practical and hygienic. And it is easy to use in a wide range of situations; i.e. patient care, private chores or on a camping holiday, when there are limited opportunities for a proper bath.



### Specifications

|                                  |   |
|----------------------------------|---|
| Base name                        | Wash Glove, Molton  |
| Brand                            | ABENA   |
| Colour                           | White   |
| Features                         | Ultrasonic sealed   |
| Single or multiple use           | Single use  |
| Material                         | Polyester   |
| Ingredient list                  | 90% Polyester, 10% Polypropylenes   |
| Length/depth                     | 23 cm   |
| Width                            | 16 cm   |
| Weight, gross                    | 0.27 kg   |
| Thickness                        | 77 g/m2   |
| Directives, regulations and acts | (EU) 2023/988   |
| Safety Instructions And Warnings | This product is for externally use only. Avoid contact with eyes.                     |
| Shelf Life                       | 60 months   |
| Shelf Life, opened               | 60 months   |
| Storage Instructions             | Store frost free in original packaging. Avoid high temperature and direct sunlight.   |
| Product Disposal Instructions    | Can be disposed of with normal household waste sorted according to local regulations. |
| Packaging Disposal Instructions  | Can be disposed of with normal household waste sorted according to local regulations. |

### Instructions for use/application

First, moisten the glove with water if you use a product from our blue range (wash with water). You can also use the glove without water, together with a product from our green range. The glove can be turned over with a slight movement, without you having to touch the unclean side. This makes the wash glove extremely hygienic. After use, remove the wash glove inside out and dispose of it together with normal household waste. For external use only.

### Packaging data

| Unit | Contains | Length  | Width   | Height  | EAN           |
|------|----------|---------|---------|---------|---------------|
| cll  | 1000 pcs | 40.5 cm | 40.5 cm | 26.5 cm | 5703538094460 |
| pck  | 50 pcs   | 26.6 cm | 20.9 cm | 6.4 cm  | 5703538094453 |
| pcs  | 1 pcs    |         |         |         |               |



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



The product is only for single-use.



The product has been tested for ingredients which can potentially cause skin irritation on people with sensitive skin. Tests conducted by dermatologists on an independent third-party laboratory.



The product has ultrasonic seals.

