



Product Main Number 70152



Pad, ABENA Let, Super Airlaid, white, with SAP

- ✓ Non-woven material with flow-through backsheet, used as a booster; 100% breathability supports healthy skin; different sizes and absorbency levels



Product Description

The ABENA Let product range is a rectangular pad solution for broad applications. All ABENA Let Airlaid Super, Normal and Maxi are used mainly as a booster pad to provide extra absorption, but can also be used as a maternity pad or functions as a secondary absorbant pad for other bodily fluids. The fully breathable pads with a non-woven backsheet help maintain healthy skin and ensure optimal comfort for the user. They are available in different sizes and absorption levels. ABENA Let Normal and Airlaid both contain superabsorbent polymer (SAP) as part of the Odour Control system.



Specifications

Base name	Pad
Brand	ABENA
Sub-Brand	Let
Size	Super Airlaid
Colour	White
Features	With SAP
Material	Nonwoven, SAP
Length/depth	32 cm
Width	12 cm
Tolerance	+/- 10%
Absorption Rothwell	140 ml
Gender	Unisex
Certificates	CE. MD. UKCA.
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745
Safety Instructions And Warnings	Keep out of reach of children
Shelf Life	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations.
Packaging Disposal Instructions	Can be recycled or incinerated.

Instructions for use/application

Rectangular pad for wide use. Mainly used as a booster together with another incontinence product, to increase overall absorbency. If it is not used as a booster,

it is recommended to wear fixation pants, for example our ABENA Fix products. Fixation pants ensure a good fit and reduce the risk of leakage. Product size recognition using colour codes and names.

Packaging data

Unit	Contains	Length	Width	Height	EAN
car	6 pck	41.5 cm	38.6 cm	27.6 cm	7331246005423
pck	100 pcs	18.4 cm	38.9 cm	14.6 cm	7331246005416



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products and medical equipment.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



Odor control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal Odor arises from the product.



The bag is 30% PCR, which means made from plastic materials used by consumers. This will allow us to reduce our plastic waste and to help take care of the environment.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



The Flow Through symbol indicates that the incontinence product can be used as a "booster", which means that the product can be placed into other incontinence products to provide extra absorbency. The booster products have a backsheet made of non-woven materials, where urine can flow through the backside of the product.



The product does not contain latex.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 9 drops.

